

TECHNICAL FEATURES


- MODEL: Fisiosonic Plus®
- COMMERCIAL CLASSIFICATION: equipment for ultrasound therapy in Class E.
- TECHNICAL CLASSIFICATION: medical device, Class I type BF.
- MEDICAL DEVICE CLASS: : IIa (Dir. 93/42/EEC, modified by the Directive 2007/47/EC).
- REFERENCES STANDARDS: IEC 60601-1, IEC 60601-2-5, IEC 60601-1-2, IEC 60601-1-6.
- SIZE: 320x430x980 mm.
- WEIGHT: ≈ 25 Kg.
- POWER: ~230V single phase (~115V available upon request).
- LINE FREQUENCY: 50-60 Hz.
- POWER CONSUMPTION: 80VA.
- LINE FUSES: 1,6A delayed.
- ULTRASOUND FREQUENCY: 1-3MHz
- POWER DENSITY ON THE TRANSDUCERS: it may be regulated continuously 0÷3 W/cm².
- OUTPUT POWER: Max. 30W (applicator Ø 55mm) - Max. 15W (applicator Ø 35mm).
- OPERATING MODE: continuous (CW) and pulsed.
 - a) Continuous mode: ultrasound energy is generated continuously and the output power density is selected by the user
 - b) Pulsed mode: modulation with adjustable frequency from 1 to 200Hz with duty-cycle ranging from 10% up to 90%.
- Visual and sound contact CONTROL SYSTEM and SELF-CALIBRATION SYSTEM.
- ENTRANCE for the connection to the electrotherapy devices (Modulo series).
- TIMER: treatment time ranging from 1 min up to 30 min., with resolution 1 min.
- ACOUSTIC SIGNAL at end of treatment.
- PRESET PROGRAMS: 40 for the main pathologies.
- FREE PROGRAMS: 100 user's programs.
- STANDARD ACCESSORIES:
 - Multifrequency applicator Ø 35mm 1/3 MHz - 5 cm² surface.
 - Multifrequency applicator Ø 55mm 1/3 MHz - 10cm² surface - attached to the articulated arm.
 - Bottle of ultrasound gel - Power supply cable - User Manual.
- OPTIONAL ACCESSORIES:
 - Multifrequency applicator Ø 21mm 1/3 MHz - 1cm² surface.
 - Rotary multidisc scanning multifrequency applicator (3 discs) 1/3 MHz - 32cm² contact surface.
 - Multidisc scanning multifrequency applicator (3 discs) 1/3 MHz - 52cm² contact surface .
 - Double multidisc scanning multifrequency applicator (6 discs) 1/3 MHz - 52cm² contact surface each applicator.

VERSION WITH GRAPHIC DISPLAY

- User Interface equipped with keyboard.
- Backlit graphic display.
- Large preset protocols database for the main pathologies.
- User programs.

TOUCH SCREEN VERSION

- Colored TFT touch screen Display 7".
- Large preset protocols database for the main pathologies.
- User programs.
- Wizard and inner Finder.
- Upgradable via USB port.

-  MARK : device complying with Directive 93/42/EEC, modified by the Directive 2007/47/EC.

 FISIOLINE® instruments are covered by UNIPOLSAI product liability insurance.



FISIOSONIC® PLUS

THE MULTIFREQUENCY ULTRASOUND IN CLASS E



MADE IN ITALY

Fisiline® srl

Borgata Molino, 29 • 12060 VERDUNO (CN) • ITALY

Tel.: +39.0172.470432-0172.470433 • Fax.: +39.0172.470891

<http://www.fisiline.com> • e-mail: fisiline@fisiline.com



Ed. 01/2020

Technical features and sizes of the devices described herein may be subject to change, based on technological evolution and constant improvements to production.

Fisiosonic Plus® is a new device for ultrasound therapy, using an innovative output stage called Class E. The mechanical efficiency obtained through the output stage in Class E is decidedly higher than the one that can be obtained by means of traditional output stages and allows the operator to issue high powers on the patient in a constant and controlled way in relation to impedance of the tissues and with a limited thermal effect.

Fisiosonic Plus® is a last generation ultrasound device that uses 1-3MHz multifrequency applicators. These frequencies can be used separately or in joint mode during the treatment.

Besides traditional methods, **Fisiosonic Plus®** has introduced two innovative functions: **sequential treatment** and **automatic frequency switching**.

The **sequential treatment** offers the operator the opportunity of using **multidisc scanning applicators** with a wide contact surface to be applied by common elastic straps. The ultrasound emission is sequential and allows the operator to carry out an automatic treatment on a wide surface with no need for a manual dedicated intervention.



Indications:

- EPICONDYLITIS
- MUSCULAR LESIONS
- IMPINGEMENT SYNDROME
- TENDINITIS
- FIBROUS ADHESIONS
- FIBROMYALGIA
- MUSCULAR TENSIONS
- MUSCLES CONTRACTURES
- ACTIVE HYPERAEMIA

The **automatic frequency switching** is obtained by the exclusive **A.F.S. (AUTOMATIC FREQUENCY SWITCHING) device**. This device makes it possible to exploit the various focalizations of 1 and 3 MHz frequencies at the same time, allowing to set different supplying times for each frequency. This feature allows the operator to obtain a focalized action.

The applicators are waterproof for underwater treatments.

The instrument can be operated in **two modes: continuous** and **pulsed**. The pulsed operating mode permits to modulate the energy (duty-cycle from 10% to 90%) in order to reduce the thermal effect and keep the peak pressure applied by the acoustic wave on the tissues steady. The output power density can be adjusted from 0 W/cm² to 3 W/cm².

Fisiosonic Plus® is equipped with **self-calibration system** and **visual and sound contact control system**. The appliance is equipped with an inlet for electrotherapy appliance connection (Modulo series); the operator can then carry out a **combined therapy** with any applicator.

Fisiosonic Plus® is equipped with a **library of preset therapeutic protocols** focused on severe and chronic pathologies. It provides the operator the **option to create customized programs** for a patient and indicated for specific pathologies or areas to be treated, for rapid and efficient therapies.

The user interface, equipped with keyboard and graphic display, has a decidedly modern and professional look.

OPTIONAL ACCESSORIES:

- A - Multifrequency applicator Ø 21mm.
- B - Rotary multidisc scanning multifrequency applicator (3 discs Ø 10 mm).
- C - Multidisc scanning multifrequency applicator (3 discs Ø 30 mm).
- D - Double multidisc scanning multifrequency applicator (6 discs Ø 30 mm).
- E - Elastic band.



VERSATILE ARTICULATED ARM FOR AUTOMATIC TREATMENTS

WATERPROOF APPLICATORS FOR UNDERWATER TREATMENTS



AVAILABLE ALSO IN TOUCH SCREEN VERSION

The interface of **Fisiosonic Plus®** in **Touch Screen** version displays the parameters on a **large color display** that allows to **remotely read** and gives to the operator the possibility to **select, modify and keep always under control** the functions and the working parameters.

A unique feature of **Fisiosonic Plus®** touch-screen version is the **Wizard**, a **step-by-step digital guide** to use the device.

Thanks to an **intuitive software**, the Wizard leads the operator in selecting the **therapeutic preset protocols** focused on **severe and chronic pathologies**.



PRESET PROTOCOLS

RAPID RESEARCH WITH THE INNER FINDER
for pathology, n. of program, body area

SELECTION FREE PROTOCOL
frequency, frequency switching, duty-cycle, duration

CREATION OF PERSONALIZED PROGRAMS
Saving data sets, 100 free storage locations, multiphase programs up to 4 phases

UPGRADABLE VIA USB PORT